

Anti-Luc7L2 Polyclonal Antibody For Wb, Elisa - Q9Y383

Item Number: CM0015559

Introduction

Rabbit polyclonal antibody against LUC7L2, validated for WB and ELISA. Targets human and mouse LUC7L2, a putative RNA-binding protein localized to nuclear speckles. Ideal for splicing and gene expression research.

[Learn More](#)

Parameter	Specification
Product Name	LUC7L2 Rabbit pAb
Remarks/Alias	CGI-59; CGI-74; LUC7B2; LUC7L2
Species	Human
GenelD (human)	51631
GenelD	51631
Immunogen	Recombinant protein
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q9Y383
Protein Weight	47kDa
Shipping	Ice bag