

Anti-Lumican (Lum) Rabbit Monoclonal Antibody For Wb, Ihc-P, Elisa - P51884

Item Number: CM0003118

Introduction

Recombinant rabbit monoclonal antibody against human Lumican (LUM), validated for WB, IHC-P, ELISA with cross-reactivity to mouse. Suitable for extracellular matrix and corneal research.

[Learn More](#)

Parameter	Value
Product Name	Lumican (LUM) Rabbit mAb
Remarks/Alias	LDC; SLRR2D; Lumican (LUM)
Species	Human
GeneID (Human)	4060
GeneID	4060
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 212-338 of human Lumican (LUM) (NP_002336.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse
SWISS	P51884
Protein Weight	38kDa
Shipping	Ice bag