

Anti-Luteinizing Hormone Subunit Beta (Lhb) Polyclonal Antibody For Wb, Ihc-P, Elisa - P01229

Item Number: CM0026537

Introduction

Rabbit polyclonal antibody targeting human LHB (luteinizing hormone beta subunit). Suitable for WB, IHC-P, ELISA; cross-reacts with mouse and rat. Ideal for reproductive endocrinology and hormone signaling studies.

[Learn More](#)

Parameter	Value
Product Name	LHB Rabbit pAb
Remarks/Alias	CGB4; HH23; LSH-B; LSH-beta; LHB
Species	Human
Gene ID (Human)	3972
Gene ID	3972
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 62-141 of human LHB (NP_000885.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P01229
Protein Weight	15kDa
Shipping	Ice bag