

## Anti-Lutropin Subunit Beta (Lhb) Polyclonal Antibody For Wb, Ihc-P, Elisa - P01229

Item Number: CM0018742

### Introduction

Rabbit polyclonal antibody against human Lutropin subunit beta (LHB) for WB, IHC-P, and ELISA. Validated in human samples; suitable for reproductive endocrinology research. Predicted MW ~15 kDa.

[Learn More](#)

| Parameter        | Value                                                                                                |
|------------------|------------------------------------------------------------------------------------------------------|
| Product Name     | LHB Rabbit pAb                                                                                       |
| Remarks/Alias    | LHB; CGB4; HH23; LSH-B; LSH-beta; lutropin subunit beta; Lutropin subunit beta                       |
| Species          | Human                                                                                                |
| GeneID (Human)   | 3972                                                                                                 |
| GeneID           | 3972                                                                                                 |
| Immunogen        | Recombinant Protein corresponding to a sequence within amino acids 21-141 of human LHB (NP_000885.1) |
| Source           | Rabbit                                                                                               |
| Category         | Polyclonal Antibodies                                                                                |
| Application      | WB, IHC-P, ELISA                                                                                     |
| Cross Reactivity | Human                                                                                                |
| SWISS            | P01229                                                                                               |
| Protein Weight   | 15kDa                                                                                                |
| Shipping         | Ice bag                                                                                              |