

Anti-Lxr Beta/Ner Rabbit Polyclonal Antibody For Wb, Elisa - P55055

Item Number: CM0026162



Introduction

KD-validated rabbit polyclonal antibody targeting human LXR beta/NER (NR1H2). Ideal for Western Blot and ELISA applications in cholesterol metabolism and transcriptional regulation studies. SWISS: P55055.

[Learn More](#)

Parameter	Value
Product Name	[KD Validated] LXR beta/NER Rabbit pAb
Remarks/Alias	NER; UNR; LXRβ; LXR-b; NER-I; RIP15; NR1H2
Species	Human
Gene ID	7376
Immunogen	Synthetic peptide corresponding to a sequence within amino acids 1-100 of human LXR beta/NER (NP_009052.4)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	P55055
Protein Weight	51kDa
Shipping	Ice bag