


Anti-Lxr α Rabbit Monoclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - Q13133

Item Number: CM0011290



Introduction

Rabbit monoclonal antibody against human LXR α (Q13133), suitable for WB, IHC-P, IF/ICC, and ELISA. Cross-reacts with mouse and rat. LXR α is a nuclear receptor key to cholesterol homeostasis regulation.

[Learn More](#)

Parameter	Value
Product Name	LXR α Rabbit mAb
Remarks/Alias	LXRA; LXR-a; RLD-1; LXR α
Species	Human
Gene ID	10062
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 1-100 of human LXR α (Q13133).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q13133
Protein Weight	50kDa
Shipping	Ice bag