


## Anti-Lxr $\alpha$ Rabbit Polyclonal Antibody For Wb And Elisa - Q13133

Item Number: CM0019495



### Introduction

High-quality rabbit polyclonal antibody targeting human oxysterols receptor LXR-alpha (Liver X Receptor alpha), validated for Western blot and ELISA. Cross-reacts with mouse and rat. Ideal for cholesterol homeostasis and nuclear receptor signaling research.

[Learn More](#)

| Parameter        | Value                                                                                                                     |
|------------------|---------------------------------------------------------------------------------------------------------------------------|
| Product Name     | LXR $\alpha$ Rabbit pAb                                                                                                   |
| Remarks/Alias    | LXRA; LXR-a; RLD-1; LXR $\alpha$                                                                                          |
| Species          | Human                                                                                                                     |
| Gene ID          | 10062                                                                                                                     |
| Immunogen        | Recombinant fusion protein containing a sequence corresponding to amino acids 160-230 of human LXR $\alpha$ (NP_005684.2) |
| Source           | Rabbit                                                                                                                    |
| Category         | Polyclonal Antibodies                                                                                                     |
| Application      | WB, ELISA                                                                                                                 |
| Cross Reactivity | Human, Mouse, Rat                                                                                                         |
| SWISS            | Q13133                                                                                                                    |
| Protein Weight   | 50kDa                                                                                                                     |
| Shipping         | Ice bag                                                                                                                   |