

Anti-Ly-6G Monoclonal Antibody For Fc - P35461

Item Number: CM0012463

Introduction

Rabbit monoclonal antibody targeting mouse Ly-6G (Gr-1), validated for flow cytometry. Ideal for immunophenotyping of granulocytes in mouse bone marrow and other tissues.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Mouse Ly-6G mAb
Remarks/Alias	Gr1; Gr-1; Ly-6G
Species	Mouse
GeneID (Human)	80740
GeneID	546644
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 27-134 of mouse Ly-6G (NP_001297367.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Mouse
SWISS	P35461
Protein Weight	28kDa
Shipping	Ice bag