

Anti-Ly-6G Rabbit Polyclonal Antibody For Wb, If-P, Elisa - P35461

Item Number: CM0023721

Introduction

Ly-6G Rabbit pAb is a polyclonal antibody validated for WB, IF-P, and ELISA. It specifically detects mouse Ly-6G, a GPI-anchored protein expressed in bone marrow, ideal for myeloid cell studies.

[Learn More](#)

Parameter	Specification
Product Name	Ly-6G Rabbit pAb
Remarks/Alias	Gr1; Gr-1; Ly-6G
Species	Mouse
GeneID	546644
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 27-105 of mouse Ly-6G (NP_001297367.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF-P, ELISA
Cross Reactivity	Mouse
SWISS	P35461
Protein Weight	14kDa
Shipping	Ice bag