

Anti-Ly6E Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - Q16553

Item Number: CM0018602

Introduction

Rabbit polyclonal antibody against human LY6E (RIG-E/SCA-2); validated for WB, IF/ICC, ELISA; cross-reacts with mouse and rat. Ideal for T-cell signaling and coronavirus entry research.

[Learn More](#)

Parameter	Value
Product Name	LY6E Rabbit pAb
Remarks/Alias	RIGE; SCA2; RIG-E; SCA-2; TSA-1; LY6E
Species	Human
GeneID (human)	4061
GeneID	4061
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 21-101 of human LY6E (NP_002337.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Mouse, Rat
SWISS	Q16553
Protein Weight	14kDa
Shipping	Ice bag