


Anti-Ly75 (Dec-205/Cd205) Rabbit Polyclonal Antibody For Wb/Elisa - O60449

Item Number: CM0017031



Introduction

Rabbit polyclonal antibody targeting human LY75 (DEC-205, CD205), raised against recombinant protein (aa 216-501). Validated for WB and ELISA with cross-reactivity to human and mouse. Suitable for endocytic receptor studies.

[Learn More](#)

Parameter	Value
Product Name	LY75 Rabbit pAb
Remarks/Alias	CD205; LY-75; CLEC13B; DEC-205; GP200-MR6; LY75
Species	Human
GeneID (Human)	4065
GeneID	4065
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 216-501 of human LY75 (NP_002340.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	O60449
Protein Weight	198kDa
Shipping	Ice bag