

Anti-Ly86 Polyclonal Antibody For Wb, Ihc-P, Elisa - O95711

Item Number: CM0030217



Introduction

CamelBio offers anti-LY86 (MD-1) rabbit polyclonal antibody validated for WB, IHC-P, and ELISA. Reacts with human and mouse LY86. Ideal for innate immunity and LPS signaling research.

[Learn More](#)

Parameter	Value
Product Name	LY86 Rabbit pAb
Remarks/Alias	MD1; MD-1; MMD-1; dJ80N2.1; LY86
Species	Human
GeneID (Human)	9450
GeneID	9450
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 21-162 of human LY86 (NP_004262.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse
SWISS	O95711
Protein Weight	18kDa
Shipping	Ice bag