


Anti-Ly96 Rabbit Polyclonal Antibody For Wb, Elisa - Q9Y6Y9

Item Number: CM0024430



Introduction

Rabbit polyclonal antibody against human LY96 (MD-2), validated for WB and ELISA; cross-reacts with mouse. An essential tool for studying LPS/TLR4 signaling and innate immunity.

[Learn More](#)

Parameter	Value
Product Name	LY96 Rabbit pAb
Remarks/Alias	MD2; MD-2; ly-96; ESOP-1; LY96
Species	Human
GeneID (Human)	23643
GeneID	23643
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 19-160 of human LY96 (NP_056179.3).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Mouse
SWISS	Q9Y6Y9
Protein Weight	19kDa
Shipping	Ice bag