


Anti-Lyn Polyclonal Antibody For Wb, Elisa - P07948

Item Number: CM0023162



Introduction

Rabbit polyclonal antibody targeting human Lyn (JTK8, p53Lyn, p56Lyn) validated for WB and ELISA; reacts with human, mouse, and rat samples. Ideal for studying B-cell signaling and immune regulation.

[Learn More](#)

Parameter	Specification
Product Name	Lyn Rabbit pAb
Remarks/Alias	JTK8; p53Lyn; p56Lyn; Lyn
Species	Human
GeneID (Human)	4067
GeneID	4067
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-240 of human Lyn (NP_002341.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P07948
Protein Weight	59kDa
Shipping	Ice bag