

Anti-Lypla1 Rabbit Polyclonal Antibody For Wb, Ihc-P, Elisa - O75608

Item Number: CM0031032



Introduction

LYPLA1 Rabbit Polyclonal Antibody validated for WB, IHC-P, and ELISA. Detects human and mouse acyl-protein thioesterase 1 (25 kDa), a key regulator of protein palmitoylation and lysophospholipid metabolism.

[Learn More](#)

Parameter	Value
Product Name	LYPLA1 Rabbit pAb
Remarks/Alias	APT1; LPL1; APT-1; LPL-I; hAPT1; LYPLA1
Species	Human
Gene ID (Human)	10434
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-230 of human LYPLA1 (NP_006321.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse
SWISS	O75608
Protein Weight	25kDa
Shipping	Ice bag