

Anti-Lyrm1 Rabbit Polyclonal Antibody For Wb And Elisa - O43325

Item Number: CM0018186

Introduction

Rabbit polyclonal antibody targeting human LYRM1, suitable for Western blot and ELISA with mouse cross-reactivity. Investigate LYRM1's role in preadipocyte proliferation and apoptosis in adipose tissue.

[Learn More](#)

Parameter	Value
Product Name	LYRM1 Rabbit pAb
Remarks/Alias	A211C6.1; LYRM1
Species	Human
GeneID (Human)	57149
GeneID	57149
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-122 of human LYRM1 (NP_065157.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Mouse
SWISS	O43325
Protein Weight	14kDa
Shipping	Ice bag