

Anti-Lysozyme (Lyz) Monoclonal Antibody For Wb, Ihc-P, Elisa - P61626

Item Number: CM0003032



Introduction

Recombinant monoclonal antibody against human Lysozyme C (LYZ) validated for WB, IHC-P, ELISA. Suitable for studies on innate immunity, antimicrobial defense, and amyloidosis research.

[Learn More](#)

Parameter	Value
Product Name	Lysozyme (LYZ) Rabbit mAb
Remarks/Alias	LZM; LYZF1; Lysozyme (LYZ)
Species	Human
GeneID (Human)	4069
GeneID	4069
Immunogen	Synthetic peptide corresponding to a sequence within amino acids 1-100 of human Lysozyme (LYZ) (NP_000230.1).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P61626
Protein Weight	17kDa
Shipping	Ice bag