

Anti-Lysozyme (Lyz) Rabbit Monoclonal Antibody For Wb, If-P, Elisa - P61626

Item Number: CM0002221



Introduction

Recombinant rabbit monoclonal antibody against human Lysozyme C (LYZ). Validated for WB, IF-P, ELISA; cross-reacts with mouse. Targets bacteriolytic enzyme in innate immunity research.

[Learn More](#)

Parameter	Value
Product Name	Lysozyme (LYZ) Rabbit mAb
Remarks/Alias	LZM; LYZF1; Lysozyme (LYZ)
Species	Human
Gene ID (Human)	4069
Immunogen	Synthetic Peptide A synthetic peptide corresponding to a sequence within amino acids 1-100 of human Lysozyme (LYZ) (P61626).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF-P, ELISA
Cross Reactivity	Human, Mouse
SWISS	P61626
Protein Weight	17kDa
Shipping	Ice bag