


Anti-Lysrs/Kars Polyclonal Antibody For Wb, If/Icc, Elisa - Q15046

Item Number: CM0017049



Introduction

Rabbit polyclonal antibody against human LysRS/KARS (KARS1), suitable for WB, IF/ICC, ELISA in human, mouse, rat samples. LysRS catalyzes tRNA aminoacylation and acts as a signaling molecule. Validated for research use.

[Learn More](#)

Parameter	Value
Product Name	LysRS/KARS Rabbit pAb
Remarks/Alias	KRS; KARS; KARS2; LEPID; CMTRIB; DEAPLE; DFN89; LysRS/KARS
Species	Human
Gene ID (Human)	3735
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 298-597 of human LysRS/LysRS/KARS (NP_005539.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q15046
Protein Weight	68kDa
Shipping	Ice bag