


Anti-Lyve1 Rabbit Monoclonal Antibody For Wb, Elisa - Q9Y5Y7

Item Number: CM0008738



Introduction

High-quality rabbit monoclonal antibody targeting human LYVE1 (lymphatic vessel endothelial hyaluronic acid receptor 1). Validated for WB and ELISA; cross-reacts with mouse and rat. Suitable for studying lymphatic endothelial cell biology and hyaluronan transport.

[Learn More](#)

Specification	Value
Product Name	LYVE1 Rabbit mAb
Remarks/Alias	HAR; XLKD1; LYVE-1; CRSBP-1; LYVE1
Species	Human
Gene ID (Human)	10894
Gene ID	10894
Immunogen	Synthetic Peptide
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9Y5Y7
Protein Weight	35kDa
Shipping	Ice bag