

Anti-Lyve1 Rabbit Monoclonal Antibody For Wb, If-P, Ihc-P, Elisa - Q8Bhc0

Item Number: CM0005050



Introduction

Mouse LYVE1-specific rabbit monoclonal antibody for WB, IF-P, IHC-P, and ELISA. Detects mouse and rat LYVE1 (34kDa), a lymphatic vessel endothelial hyaluronic acid receptor involved in HA transport and cell adhesion.

[Learn More](#)

Parameter	Value
Product Name	LYVE1 Rabbit mAb
Remarks/Alias	Xlkd1; Lyve-1; Crsbp-1; 1200012G08Rik; Lyve1
Species	Mouse
GeneID (human)	10894
GeneID	114332
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 30-228 of mouse Lyve1. (NP_444477.2).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF-P, IHC-P, ELISA
Cross Reactivity	Mouse, Rat
SWISS	Q8BHC0
Protein Weight	34kDa
Shipping	Ice bag