

Anti-Lyve1 Rabbit Polyclonal Antibody For Wb And Elisa - Q8Bhc0

Item Number: CM0018907



Introduction

Lyve1 Rabbit Polyclonal Antibody is a rabbit pAb raised against mouse Lyve1 (aa 30-228), validated for Western blot and ELISA, and cross-reacts with mouse and rat. Ideal for studying lymphatic endothelial hyaluronan receptor 1 function and hyaluronan transport.

[Learn More](#)

Parameter	Value
Product Name	Lyve1 Rabbit pAb
Remarks/Alias	Xlkd1; Lyve-1; Crsbp-1; 1200012G08Rik; Lyve1
Species	Mouse
GeneID	114332
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 30-228 of mouse Lyve1 (NP_444477.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Mouse, Rat
SWISS	Q8BHC0
Protein Weight	34kDa
Shipping	Ice bag