

Anti-Lyxl6 Polyclonal Antibody For Wb And Elisa - 075951

Item Number: CM0016900

Introduction

Rabbit polyclonal antibody targeting human LYXL6, validated for Western blot (WB) and ELISA. Suitable for fertility and antimicrobial research, detecting a ~17 kDa protein from testis and sperm samples.

[Learn More](#)

Parameter	Value
Product Name	LYXL6 Rabbit pAb
Remarks/Alias	LYC1; LYXB; UNQ754; PRO1485; TKAL754; HEL-S-6a; LYXL6
Species	Human
Gene ID (Human)	57151
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 20-148 of human LYXL6 (NP_065159.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	O75951
Protein Weight	17kDa
Shipping	Ice bag