

Anti-Lztr1 Rabbit Polyclonal Antibody For Wb, Elisa - Q8N653

Item Number: CM0024950



Introduction

LZTR1 rabbit polyclonal antibody validated for WB and ELISA. Detects human, mouse, and rat LZTR1, a substrate-specific adapter of the BCR E3 ubiquitin ligase complex regulating Ras-MAPK signaling. Ideal for oncology and neurodevelopmental disease research.

[Learn More](#)

Parameter	Value
Product Name	LZTR1 Rabbit pAb
Remarks/Alias	NS2; NS10; BTBD29; LZTR-1; SWNTS2; LZTR1
Species	Human
GeneID (human)	8216
GeneID	8216
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 300-545 of human LZTR1 (NP_006758.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q8N653
Protein Weight	95kDa
Shipping	Ice bag