

Anti-Lzts2 Rabbit Polyclonal Antibody For Wb And Elisa - Q9Brk4

Item Number: CM0014863

Introduction

Rabbit polyclonal antibody targeting human LZTS2 (LAPSER1), a negative regulator of microtubule severing and Wnt signaling. Validated for WB and ELISA, cross-reacts with mouse and rat. Suitable for cell cycle and cancer research.

[Learn More](#)

Parameter	Value
Product Name	LZTS2 Rabbit pAb
Remarks/Alias	LAPSER1; LZTS2
Species	Human
GeneID	84445
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 470-669 of human LZTS2 (NP_115805.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9BRK4
Protein Weight	73kDa
Shipping	Ice bag