


Anti-MacroH2A.1 Polyclonal Antibody For Wb, If/Icc, Elisa - 075367

Item Number: CM0021307



Introduction

Rabbit polyclonal antibody targeting macroH2A.1 (H2AFY), a variant histone involved in transcriptional repression, X inactivation, and NAD⁺ metabolism. Validated for WB, IF/ICC, ELISA in human, mouse, rat.

[Learn More](#)

Parameter	Value
Product Name	[KO Validated] macroH2A.1 Rabbit pAb
Remarks/Alias	H2A.y; H2A/y; H2AFY; mH2A1; H2AF12M; MACROH2A1.1; macroH2A1.2; .1
Species	Human
Gene ID (Human)	9555
Gene ID	9555
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 123-372 of human macroH2A.1 (NP_613258.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	075367
Protein Weight	39kDa
Shipping	Ice bag