


Anti-MacroH2A.1 Rabbit Polyclonal Antibody For Wb, Elisa - O75367

Item Number: CM0021345



Introduction

Rabbit polyclonal antibody targeting human macroH2A.1 (H2AFY), validated for WB and ELISA. Recognizes human and rat isoforms; core histone variant involved in transcription repression and chromatin remodeling.

[Learn More](#)

Parameter	Value
Product Name	macroH2A.1 Rabbit pAb
Remarks/Alias	H2A.y; H2A/y; H2AFY; mH2A1; H2AF12M; MACROH2A1.1; macroH2A1.2; macroH2A.1
Species	Human
GeneID	9555
Immunogen	Synthetic peptide A synthetic peptide corresponding to a sequence within amino acids 200-300 of human macroH2A.1 (NP_613258.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Rat
SWISS	O75367
Protein Weight	39kDa
Shipping	Ice bag