


Anti-Mad1L1 Polyclonal Antibody For Wb, Elisa - Q9Y6D9

Item Number: CM0020113



Introduction

Rabbit polyclonal antibody targeting human MAD1L1 for WB and ELISA; cross-reacts with rat. Suitable for spindle assembly checkpoint, chromosome segregation, and mitotic checkpoint signaling research.

[Learn More](#)

Parameter	Value
Product Name	MAD1L1 Rabbit pAb
Remarks/Alias	MAD1; MVA7; PIG9; TP53I9; TXBP181; MAD1L1
Species	Human
GeneID (Human)	8379
GeneID	8379
Immunogen	Synthetic peptide corresponding to a sequence within amino acids 1-100 of human MAD1L1 (NP_003541.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Rat
SWISS	Q9Y6D9
Protein Weight	83kDa
Shipping	Ice bag