


Anti-Mad1/Mad1L1 Rabbit Monoclonal Antibody For Wb, Elisa - Q9Y6D9

Item Number: CM0010956



Introduction

Rabbit monoclonal antibody targeting human MAD1/MAD1L1 (MAD1, TXBP181), validated for Western blot and ELISA. Detects both human and mouse MAD1; ideal for spindle checkpoint and mitosis research.

[Learn More](#)

Parameter	Value
Product Name	MAD1/MAD1L1 Rabbit mAb
Remarks/Alias	MAD1; MVA7; PIG9; TP53I9; TXBP181; MAD1/MAD1L1
Species	Human
GeneID (Human)	8379
GeneID	8379
Immunogen	Synthetic Peptide
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q9Y6D9
Protein Weight	83kDa
Shipping	Ice bag