


Anti-Mad2B/Mad2L2 Polyclonal Antibody For Wb, If/Icc, Elisa - Q9UI95

Item Number: CM0013781



Introduction

Rabbit polyclonal antibody against human MAD2B/MAD2L2 (REV7), validated for WB, IF/ICC, ELISA. Cross-reacts with human, mouse, rat. Useful for DNA damage response and cell cycle studies.

[Learn More](#)

Parameter	Value
Product Name	MAD2B/MAD2L2 Rabbit pAb
Remarks/Alias	REV7; FANCV; MAD2B; POLZ2; MAD2B/MAD2L2
Species	Human
GeneID (Human)	10459
GeneID	10459
Immunogen	Recombinant protein
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9UI95
Protein Weight	24kDa
Shipping	Ice bag