


Anti-Mad2B/Mad2L2 Polyclonal Antibody For Wb - Q9Ui95

Item Number: CM0026857



Introduction

Rabbit polyclonal antibody against MAD2B/MAD2L2 (REV7), validated for WB and ELISA in human and mouse. This mitotic spindle checkpoint protein is involved in DNA damage repair and translesion synthesis.

[Learn More](#)

Parameter	Value
Product Name	MAD2B/MAD2L2 Rabbit pAb
Remarks/Alias	REV7; FANCV; MAD2B; POLZ2; MAD2B/MAD2L2
Species	Human
GeneID (Human)	10459
GeneID	10459
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-211 of human MAD2B/MAD2L2 (NP_006332.3)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q9UI95
Protein Weight	24kDa
Shipping	Ice bag