


Anti-Mad2L1Bp Rabbit Monoclonal Antibody For Ihc-P, If/Icc, Elisa - Q15013

Item Number: CM0002860



Introduction

Recombinant rabbit monoclonal antibody against human MAD2L1BP (p31(comet)), validated in IHC-P, IF/ICC, and ELISA. Reacts with human and mouse. Suitable for spindle checkpoint and cell cycle research.

[Learn More](#)

Parameter	Value
Product Name	MAD2L1BP Rabbit mAb
Remarks/Alias	CMT2; MAD2L1BP
Species	Human
Gene ID (Human)	9587
Immunogen	Synthetic Peptide
Source	Rabbit
Category	Monoclonal Antibodies
Application	IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q15013
Protein Weight	31kDa
Shipping	Ice bag