

Anti-Mad2/Mad2L1 Polyclonal Antibody For Ihc-P, If/Icc, Elisa - Q13257

Item Number: CM0019242



Introduction

Rabbit polyclonal antibody targeting human MAD2/MAD2L1, a critical mitotic checkpoint protein. Validated in IHC-P, IF/ICC, and ELISA with cross-reactivity to human, mouse, and rat. Suitable for cell cycle and cancer research.

[Learn More](#)

Parameter	Value
Product Name	MAD2/MAD2L1 Rabbit pAb
Remarks/Alias	MAD2; HSMAD2; MAD2/MAD2L1
Species	Human
GeneID (Human)	4085
GeneID	4085
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-205 of human MAD2/MAD2/MAD2L1 (NP_002349.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q13257
Protein Weight	24kDa
Shipping	Ice bag