

Anti-Mad2/Mad2L1 Rabbit Polyclonal Antibody For Wb, Elisa - Q13257

Item Number: CM0023365



Introduction

Rabbit polyclonal antibody against human MAD2/MAD2L1, a mitotic spindle checkpoint protein. Validated for WB and ELISA, cross-reacts with mouse and rat. Ideal for cell cycle and mitosis research.

[Learn More](#)

Parameter	Value
Product Name	MAD2/MAD2L1 Rabbit pAb
Remarks/Alias	MAD2; HSMAD2; MAD2/MAD2L1
Species	Human
GeneID (Human)	4085
GeneID	4085
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-205 of human MAD2/MAD2L1 (NP_002349.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q13257
Protein Weight	24kDa
Shipping	Ice bag