

Anti-Mad3/Mxd3 Rabbit Monoclonal Antibody For Wb, Elisa - Q9Bw11

Item Number: CM0002718

Introduction

Rabbit monoclonal antibody against human MAD3/MXD3, validated for WB and ELISA. Recognizes a 23kDa transcriptional repressor that antagonizes MYC activity, suitable for cancer and transcription research.

[Learn More](#)

Parameter	Value
Product Name	MAD3/MXD3 Rabbit mAb
Remarks/Alias	MYX; MAD3; BHLHC13; MAD3/MXD3
Species	Human
GeneID (Human)	83463
GeneID	83463
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 100-190 of human MAD3/MXD3 (Q9BW11)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	Q9BW11
Protein Weight	23kDa
Shipping	Ice bag