

Anti-Maf1 Rabbit Polyclonal Antibody For Wb, Elisa - Q9H063

Item Number: CM0015668



Introduction

High-quality Anti-MAF1 Rabbit Polyclonal Antibody validated for WB and ELISA. Targets human and mouse MAF1, a repressor of RNA polymerase III transcription. Predicted molecular weight 29kDa.

[Learn More](#)

Parameter	Value
Product Name	MAF1 Rabbit pAb
Remarks/Alias	MAF1; MAF1 homolog
Species	Human
Gene ID	84232
Immunogen	Recombinant protein corresponding to amino acids 1-256 of human MAF1 (NP_115648.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q9H063
Protein Weight	29kDa
Shipping	Ice bag