


## Anti-Mafa Polyclonal Antibody For Wb And Elisa - Q8Nhw3

Item Number: CM0022342



### Introduction

Anti-MAFA rabbit polyclonal antibody validated for WB, ELISA; cross-reacts with mouse. Targets transcription factor MafA (37 kDa) involved in insulin gene expression. Suitable for diabetes research.

[Learn More](#)

| Parameter        | Value                                                                                                             |
|------------------|-------------------------------------------------------------------------------------------------------------------|
| Product Name     | MAFA Rabbit pAb                                                                                                   |
| Remarks/Alias    | INSDM; hMafA; RIPE3b1; MAFA                                                                                       |
| Species          | Human                                                                                                             |
| GeneID (human)   | 389692                                                                                                            |
| GeneID           | 389692                                                                                                            |
| Immunogen        | Recombinant fusion protein containing a sequence corresponding to amino acids 269-353 of human MAFA (NP_963883.2) |
| Source           | Rabbit                                                                                                            |
| Category         | Polyclonal Antibodies                                                                                             |
| Application      | WB, ELISA                                                                                                         |
| Cross Reactivity | Mouse                                                                                                             |
| SWISS            | Q8NHW3                                                                                                            |
| Protein Weight   | 37kDa                                                                                                             |
| Shipping         | Ice bag                                                                                                           |