

Anti-Maff Rabbit Polyclonal Antibody For Wb And Elisa - Q9ULx9

Item Number: CM0017895

Introduction

This anti-MAFF rabbit polyclonal antibody detects human transcription factor MafF, validated for Western blotting and ELISA. Suitable for research on gene regulation, oxytocin receptor expression, and cellular stress responses.

[Learn More](#)

Parameter	Value
Product Name	MAFF Rabbit pAb
Remarks/Alias	U-MAF; hMafF; MAFF
Species	Human
Gene ID	23764
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 75-164 of human MAFF (NP_036455.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	Q9ULX9
Protein Weight	18kDa
Shipping	Ice bag