

Anti-Mafk Rabbit Polyclonal Antibody For Wb, Elisa - O60675

Item Number: CM0022641

Introduction

Rabbit polyclonal antibody targeting human transcription factor MafK (NF-E2 p18), validated for WB and ELISA. Useful for studying transcriptional regulation and erythroid differentiation.

[Learn More](#)

Parameter	Value
Product Name	MAFK Rabbit pAb
Remarks/Alias	P18; NFE2U; MAFK
Species	Human
Gene ID	7975
Immunogen	Synthetic peptide within amino acids 1-50 of human MAFK (NP_002351.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	O60675
Protein Weight	18kDa
Shipping	Ice bag