

Anti-Mag Monoclonal Antibody For Wb, Elisa - P20916

Item Number: CM0010473

Introduction

Rabbit monoclonal antibody targeting myelin-associated glycoprotein (MAG/Siglec-4a). Validated for Western blot and ELISA. Cross-reacts with human, mouse, and rat. Ideal for neuroscience and myelination research.

[Learn More](#)

Parameter	Value
Product Name	MAG Rabbit mAb
Remarks/Alias	GMA; S-MAG; SPG75; SIGLEC4A; SIGLEC-4A; MAG
Species	Human
GeneID (Human)	4099
GeneID	4099
Immunogen	Synthetic Peptide A synthetic peptide corresponding to a sequence within amino acids 1-100 of human MAG (P20916).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P20916
Protein Weight	69kDa
Shipping	Ice bag