

Anti-Mag Monoclonal Antibody For Wb, Ihc-P, If-F, Ip, Elisa - P20916

Item Number: CM0006245

Introduction

High-quality rabbit monoclonal antibody against human MAG (myelin-associated glycoprotein). Validated for WB, IHC-P, IF-F, IP, ELISA; cross-reacts with human, mouse, rat. Ideal for neurobiology and myelination research.

[Learn More](#)

Parameter	Value
Product Name	MAG Rabbit mAb
Remarks/Alias	GMA; S-MAG; SPG75; SIGLEC4A; SIGLEC-4A
Species	Human
GenelD (Human)	4099
GenelD	4099
Immunogen	Recombinant protein
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF-F, IHC-P, IP, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P20916
Protein Weight	69kDa
Shipping	Ice bag