

Anti-Mag (Myelin-Associated Glycoprotein) Polyclonal Antibody For Wb And Elisa - P20916

Item Number: CM0025248

Introduction

Rabbit polyclonal antibody targeting human MAG (myelin-associated glycoprotein, Siglec-4a). Validated for WB and ELISA with cross-reactivity to mouse and rat. Ideal for neuroscience and myelination research.

[Learn More](#)

Parameter	Value
Product Name	MAG Rabbit pAb
Remarks/Alias	GMA; S-MAG; SPG75; SIGLEC4A; SIGLEC-4A; MAG
Species	Human
GeneID (Human)	4099
GeneID	4099
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 20-220 of human MAG (P20916)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P20916
Protein Weight	69kDa
Shipping	Ice bag