


Anti-Mage-1/Mage1A Monoclonal Antibody For Wb And Elisa - P43355

Item Number: CM0003661



Introduction

Rabbit monoclonal antibody specific to MAGE-1/MAGE1A, validated for Western blotting and ELISA. Recognizes human and mouse MAGE-1, a cancer/testis antigen involved in transcriptional regulation and tumor immunity.

[Learn More](#)

Parameter	Value
Product Name	MAGE-1/MAGE1A Rabbit mAb
Remarks/Alias	CT1.1; MAGE1; MAGE-1/MAGE1A
Species	Human
GeneID (Human)	4100
GeneID	4100
Immunogen	Synthetic Peptide; A synthetic peptide corresponding to a sequence within amino acids 200-309 of human MAGE-1/MAGE1A (P43355).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	P43355
Protein Weight	34kDa
Shipping	Ice bag