


Anti-Mage-1/Mage1A Rabbit Polyclonal Antibody For Wb - P43355

Item Number: CM0030314



Introduction

Rabbit polyclonal antibody against human MAGE-1/MAGE1A (MAGEA1) for WB, IHC-P, IF/ICC, and ELISA. Recognizes a cancer/testis antigen expressed in melanoma, lung, and breast carcinoma but not normal tissues except testis.

[Learn More](#)

Parameter	Value
Product Name	MAGE-1/MAGE1A Rabbit pAb
Remarks/Alias	CT1.1; MAGE1; MAGE-1/MAGE1A
Species	Human
GeneID (Human)	4100
GeneID	4100
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-309 of human MAGE-1/MAGE1A (NP_004979.3)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human
SWISS	P43355
Protein Weight	34kDa
Shipping	Ice bag