


Anti-Magea11 Rabbit Polyclonal Antibody For Wb, Elisa - P43364

Item Number: CM0015779



Introduction

MAGEA11 Rabbit Polyclonal Antibody for WB, ELISA. Detects human MAGEA11 (melanoma-associated antigen 11), an androgen receptor coregulator expressed in cancers and normal tissues.

[Learn More](#)

Parameter	Value
Product Name	MAGEA11 Rabbit pAb
Remarks/Alias	CT1.11, MAGE11, MAGE-11, MAGEA-11, MAGEA11
Species	Human
Gene ID (Human)	4110
Gene ID	4110
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-110 of human MAGEA11 (NP_005357.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	P43364
Protein Weight	48kDa
Shipping	Ice bag