


Anti-Magea4 Rabbit Polyclonal Antibody For Wb, Elisa - P43358

Item Number: CM0016277



Introduction

MAGEA4 Rabbit pAb (P43358) for WB & ELISA, cross-reactive with human, mouse, and rat. Recognizes 35 kDa melanoma-associated antigen 4. For cancer research and cell cycle studies.

[Learn More](#)

Parameter	Value
Product Name	MAGEA4 Rabbit pAb
Remarks/Alias	CT1.4; MAGE4; MAGE4A; MAGE4B; MAGE-41; MAGE-X2; MAGEA4
Species	Human
GeneID (human)	4103
GeneID	4103
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 28-317 of human MAGEA4 (NP_002353.3)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P43358
Protein Weight	35kDa
Shipping	Ice bag