

Anti-Magea6 Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - P43360

Item Number: CM0028583

Introduction

MAGEA6 Rabbit pAb is a polyclonal antibody validated for WB, IF/ICC, and ELISA, targeting human and mouse melanoma-associated antigen 6 (CT1.6). Ideal for studying ubiquitin ligase regulation and tumor progression, supplied as liquid.

[Learn More](#)

Parameter	Value
Product Name	MAGEA6 Rabbit pAb
Remarks/Alias	CT1.6, MAGE6, MAGE3B, MAGE-3b, MAGEA6
Species	Human
Gene ID (Human)	4105
Immunogen	Recombinant protein
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse
SWISS	P43360
Protein Weight	35kDa
Shipping	Ice bag