

Anti-Maged1/Nrage Rabbit Monoclonal Antibody For Wb, Elisa - Q9Y5V3

Item Number: CM0006697

Introduction

Rabbit monoclonal antibody targeting human MAGED1/NRAGE (Q9Y5V3). Validated for WB and ELISA; reacts with human MAGED1/NRAGE. Ideal for apoptosis, circadian rhythm, and cancer research.

[Learn More](#)

Parameter	Value
Product Name	MAGED1/NRAGE Rabbit mAb
Remarks/Alias	NRAGE; DLXIN-1; MAGED1/NRAGE
Species	Human
GeneID (Human)	9500
GeneID	9500
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 640-778 of human MAGED1/NRAGE. (NP_008917.3)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	Q9Y5V3
Protein Weight	86 kDa
Shipping	Ice bag