

Anti-Maged1/Nrage Rabbit Polyclonal Antibody For Wb, Ihc-P, Elisa - Q9Y5V3

Item Number: CM0015269

Introduction

Rabbit polyclonal antibody targeting MAGED1/NRAGE, validated for WB, IHC-P, ELISA in human, mouse, rat. Suitable for apoptosis and circadian rhythm research studies.

[Learn More](#)

Parameter	Value
Product Name	MAGED1/NRAGE Rabbit pAb
Remarks/Alias	NRAGE; DLXIN-1; MAGED1/NRAGE
Species	Human
Gene ID (Human)	9500
Gene ID	9500
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 635-834 of human MAGED1/NRAGE (NP_001005333.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9Y5V3
Protein Weight	86kDa
Shipping	Ice bag