

Anti-Mageh1 Polyclonal Antibody For Wb, Elisa - Q9H213

Item Number: CM0023758

Introduction

MAGEH1 Rabbit pAb targets melanoma-associated antigen H1 (APR-1, MAGE-H1, Restin), a tumor antigen and phosphoprotein. Suitable for WB and ELISA, with mouse cross-reactivity. Ideal for apoptosis and cancer research.

[Learn More](#)

Parameter	Value
Product Name	MAGEH1 Rabbit pAb
Remarks/Alias	APR1; APR-1; MAGEH; MAGEH1
Species	Human
GeneID	28986
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 10-90 of human MAGEH1 (NP_054780.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Mouse
SWISS	Q9H213
Protein Weight	24kDa
Shipping	Ice bag